

ALKALINE PHOSPHATASE KINETIC-DGKC METHOD



REF DESCRIPTION

CZ901U ALKALINE PHOSPHATASE

PACKAGE SIZES

CZ901U-1	50 ml	CZ901U-4	10 X 100 ml
CZ901U-2	100 ml	CZ901U-5	2 X 500 ml
CZ901U-3	10 X 20 ml	CZ901U-6	5 L

INTENDED USE

This reagent is intended for *in vitro* quantitative determination of Alkaline Phosphatase (ALP) in serum or plasma.

INTENDED USER: Professional Use Only.

METHOD

KINETIC-DGKC METHOD

CLINICAL SIGNIFICANCE

Alkaline Phosphatase (ALP) is widely distributed throughout the body, but clinically important one for diagnostic reasons is in bone, liver, placenta & intestine. Growing bone is associated with the release of ALP and so in childhood the level of ALP is around 3 times of that of adult. During pregnancy in 2nd & 3rd trimester the enzyme rises considerably due to placenta releasing ALP. It can be used to examine placental function. Elevated levels are seen in bone diseases, e.g. Paget's disease, rickets,

osteoplastic metastatic & in obstructive disease of biliary tract.

Decreased levels are rarely seen. E.g. in Vitamin A resistant rickets.

PRINCIPLE

Alkaline Phosphatase (ALP) catalyses the hydrolysis of P-Nitrophenyl phosphate (p-Npp) at pH 9.8, liberating p-Nitrophenol and phosphate, according to the following reaction:

ALP
p-Npp +H2O -----→ P-Nitrophenol + Phosphate

The rate of p-Nitrophenol formation, measured photometrically, is proportional to the catalytic concentration of alkaline Phosphatase present in the sample

REAGENT COMPOSITION

ALP BUFFER R1

Diethanolamine 1.0 mol/l Magnesium chloride 0.5 mmol/l

ALP SUBSTRATE R2

P-Nitrophenylphosphate 10 mmol/l

REAGENT PREPARATION SUBSTRATE START

R1 and R2 are ready-to-use and stable until expiry date if contamination is avoided and stored at 2-8 $^{\circ}\text{C}$

SAMPLE START

Mix 4 parts of R1 + 1 Part of R2 = Mono reagent Stability of mono reagent: 4 Weeks at 2-8°C, 5 days at 15-25°C

SPECIMEN

Serum, heparinized plasma

PRECAUTION

- 1- The reagents contain sodium azide as preservative. Do not swallow and avoid contact with skin and mucous membranes.
- 2- During the reaction P-Nitrophenol is produced. This is poisonous when inhaled, swallow or when absorbed through skin. If the reaction mixture comes in contact with skin or mucous membranes wash copiously with water.

3- To avoid contamination, use clean laboratory wares. Avoid direct exposure of reagent to light.

ASSAY

Wavelength : 405 nm
Cuvette : 1 cm light path
Temperature : 25°C/30°C/37°C
Measurement : Against reagent blank

PROCEDURE

SURSTRATE START

30B3TRATE START				
Reagent 1 Buffer	1000 μL			
Sample	20 μL			
Mix incubates for approx 1 min, then add				
D	250			
Reagent 2 Substrates	250 μL			
CARADIFICTART				

SAMPLE START

Mono reagent (R1+R2)	1000 μL
Sample	20 μL

READING FOR BOTH

Mix and read absorbance after 1 min and start stop watch. Read absorbance again after 1, 2, 3 min.

CALCULATION

SUBSTRATE START

ALP activity U/L = Δ A/min. X 3433

SAMPLE START

ALP activity U/L = $\Delta A/min$. X 2757

LINEARITY

The reagent is linear upto 700 U/L

If the activities exceed 700 U/L, mix 50 μL of sample with 200 μL of 0.9% NaCl solution and multiply the result by 5.

2500


NORMAL RANGE

	25 C	30 C	3/ C
Women U/L	40-190	48-223	64-306
Men U/L	50-190	60-223	80-306
Children up to 15 years	400	488	644
U/L. up to			

Each laboratory should establish reference ranges for its own patients population.

QUALITY CONTROL

All control sera with values determined by this method can be used. **SYMBOL ON LABELS**



BIBILOGRAPHY

- 1. Z Klin. Chem. Klin biochem 8 (1970) 658; 10 (1972) 182
- 2. Kubler, W., symp. D. Deutsch. Ges. Fur lab. Med. 8 (1973)
- 3. Weishaar, D., Grossau, 4., U. Fadel: B.Med.Welt 26 (1975), 387-390
- 4. Thomas, L, labor U.Duag, 1. Adfl.S.64 Die .Med verlag, Marb



CRESCENT DIAGNOSTICS FACTORY

Phase-3, Jeddah Industrial City, Po.Box.9939, Jeddah-21423, Saudi Arabia Tel:+966 12 6080425, Fax:+966 12 6080420, e-mail:info@crescentksa.com

CZ901U-ALP R2/15-12-2020

